


## MonoCap HighResolution 2000 Series

Optimized for Identification of Peptides/Proteins for Proteome Research



MonoCap HighResolution 2000 is a 2 meter length monolithic silica capillary column which is designed for identifying extremely high number of peptides/proteins for proteome research via LC-MS/MS.

GL Sciences' MonoCap capillary columns, created synthetically via sol-gel method, and an octadecyl silane chemically bonded, has a very uniform three dimensional structure that shows excellent reproducibility from batch-to-batch. The solid structure of GL Sciences' monolithic silica eliminates the need for frits or filters at the ends of the column, thereby reducing dead volume that might otherwise lead to band broadening or sample recovery. The high porosity of our monolithic silica allows high flow rates to be used without loss of resolution or creation of high operating pressure. An optimized balance of through-pores and meso-pores provides the critically important combination of efficiency, separation speed, large volume sample-loading, and small volume sample-recovery.

MonoCap HighResolution provide extremely high efficiency, delivering over 200,000 plates for a 2,000 mm length column. The MonoCap HighResolution Ultra type deliver over 300,000 plates.

### Physical Properties

Product Description	Bonded Phase	Meso-pore	End-capping	Max. Operating Pressure
MonoCap C18 HighResolution 2000	Octadecyl Groups	15 nm	Yes	35 MPa
MonoCap C18 HighResolution Ultra 2000	Octadecyl Groups	11 nm	Yes	35 MPa
MonoCap HILIC-UP HighResolution 2000	Ureidopropyl Groups	12 nm	None	35 MPa

Based on monolithic technology, Merck KGaA, Darmstadt, Germany.




## Discover New Peptides/Proteins

As proven below, MonoCap C18 HighResolution 2000 mm length column identifies simply more peptides/proteins compared to those traditional particle packed capillary HPLC columns.


**MS:** LTQ-Orbitrap XL (Mascot Search)  
**Sample:** THP-1 Cell Lysate Tryptic Digest, 5 µg

Column Name	Number of Identified Proteins in average	Analysis Time
<b>MonoCap C18 HighResolution 2000 0.1 mm I.D. x 2000 mm</b>	<b>2,087 (2013, 2116, 2131)</b>	<b>10 Hours</b>
Particle packed column 0.1 mm I.D. x 150 mm	680 (685, 679, 675)	2 Hours

### [1] MonoCap C18 HighResolution 2000 (2000 mm x 0.1 mm I.D.)



### [2] Particle Packed column (3 µm, 150 mm x 0.075 mm I.D.)



#### Conditions


System	: GLS Capillary HPLC system	Flow Rate	: [1] 0.5 µL/min
Column	: [1] MonoCap C18 High Resolution 2000 (2000 mm x 0.1 mm I.D.)		: [2] 0.3 µL/min
	: [2] Particle packed column (3 µm, 150 mm x 0.075 mm I.D.)	Injection Vol.:	: 5 µL
Trap column	: MonoCap C18 Trap Column (50 mm x 0.075 mm I.D.)	Detection	: MS (TIC m / z 500-1500)
Eluent	: A) 0.1 %HCOOH in CH <sub>3</sub> CN	Sample	: Tryptic digest of proteins
	: B) 0.1 %HCOOH in H <sub>2</sub> O		
	[1] A / B = 10 / 90 - 600 min - 45 / 55		
	[2] A / B = 10 / 90 - 180 min - 45 / 55		

● ● ● are analytes having the same molecular weight


## For Identifying Highly Hydrophilic, Hydrophobic Peptides/Proteins

MonoCap HILIC-UP is an important addition to the MonoCap C18 HighResolution 2000 column series. MonoCap HILIC-UP can retain highly hydrophilic peptides/proteins which may lead to discovering new peptides/proteins where a C18 phase couldn't identify.

In HILIC, the higher the organic concentration, the greater the retention of more polar analytes. One of the biggest benefit of HILIC mode is, a high organic solvent concentration of the mobile phase will lead to a high sensitivity LC-MS/MS analysis.



Total Ion Chromatogram



Base Peak Chromatogram

### Remarks:

#### Results of MonoCap C18 HighResolution 2000

Number of Identified Peptides: 8,358

Number of Identified Proteins: 1,992

Gradient Program: 4 hrs

#### Results of MonoCap HILIC-UP HighResolution 2000

Number of Identified Peptides: 7,194 (14,736 PSM\*)

Number of Identified Proteins: 2,201

\* Peptide Spectrum Match

### Conditions

Column : MonoCap HILIC-UP High Resolution 2000  
Eluent : A) CH<sub>3</sub>CN : H<sub>2</sub>O=10/90 (0.5% CH<sub>3</sub>COOH)  
B) CH<sub>3</sub>CN : H<sub>2</sub>O=95/5 (0.5% CH<sub>3</sub>COOH)  
A/B=0/100-(240 min)-20/80-(10 min)-100/0-(10 min)-100/0  
Flow Rate : 0.5 µL  
Injection Vol. : 1 µL (1 mg/mL)  
Detection : TIC MS (m/z 300-1500)  
Sample : Tryptic Digest of Hela Cell Lysate, 5 ug

### Reference:

Hydrophilic Interaction Chromatography Using a Meter-Scale Monolithic silica capillary Column for Proteomics LC-MS,  
K Horie et al. *Anal. Chem.* 2014, 86, 3817-3824

## References

---

1. M.H.M. van de Meent et al.  
Improvement of the liquid-chromatographic analysis of protein tryptic digests by the use of long-capillary monolithic columns with UV and MS detection, *Anal Bioanal Chem*, 2007,388, 195-200
2. Mio Iwasaki et al.  
One-Dimensional Capillary Liquid Chromatographic Separation Coupled with Tandem Mass Spectrometry Unveils the Escherichia coli Proteome on a Microarray Scale, *Anal. Chem.* 2010, 82, 2616-2620
3. Mio Iwasaki et al.  
Human Proteome analysis by using reversed phase monolithic silica capillary columns with enhanced sensitivity, *J Chromatogr A* 2012, 1228, 292-297
4. Ryota Yamana et al.  
Rapid and deep profiling of human induced pluripotent stem cell proteome by one-shot NanoLC-MS/MS analysis with meter-scale monolithic silica columns, *J Proteome Res.* 2013, 12, 214-21
5. Mari Ogawa-Ohnishi et al.  
Identification of three hydroxyproline O-arabinosyltransferases in *Arabidopsis thaliana*, *Nature Chem. Biol.* 2013, 9, 726-730
6. Satoru Okamoto et al.  
Root-derived CLE glycopeptides control nodulation by direct binding to HAR1 receptor kinase, *Nature Commun.* 2013,4, 2191
7. Kanta Horie et al.  
Hydrophilic interaction chromatography using a meter-scale monolithic silica capillary column for proteomics LC-MS, *Anal. Chem.* 2014, 86, 3817-3824

## Ordering Information

### MonoCap C18 HighResolution Ultra 2000

- End-fittings are not included.
- A column connection kit is available separately to ensure proper connections.
- Please refer to the below ordering information.

Description	I.D. (mm)	Length (mm)	Qty	Cat. No.
MonoCap C18 HighResolution Ultra 2000	10 µm	2000	1/pk	5020-10018

### MonoCap C18 HighResolution 2000

- End-fittings are not included.
- A column connection kit is available separately to ensure proper connections.
- Please refer to the below ordering information.

Description	I.D. (mm)	Length (mm)	Qty	Cat. No.
MonoCap C18 HighResolution 2000	10 µm	2000	1/pk	5020-10015

### MonoCap HILIC-UP HighResolution 2000

- End-fittings are not included.
- A column connection kit is available separately to ensure proper connections.
- Please refer to the below ordering information.

Description	I.D. (mm)	Length (mm)	Qty	Cat. No.
MonoCap HILIC-UP HighResolution 2000	10 µm	2000	1/pk	5020-10019

### Connection Kit for MonoCap HighResolution 2000

- A dedicated connection kit for MonoCap C18 High Resolution 2000.
- Use this connection kit when connecting the column directly to the system.

Description	Qty	Cat. No.
1/16" PEEK Ferrule, SUS Nut, Sleeve, 2 pcs each.	1/pk	5020-10017
1/32" PEEK Ferrule, SUS Nut, Sleeve, 2 pcs each.		



### Zero Dead Volume Union

- Connect the tubing from the system to this union and install the column to achieve zero dead volume.

Description	Orifice Size	Qty	Cat. No.
U-435	0.25 mm	1/pk	6010-72352
U-411	178 µm	1/pk	6010-72351




# MonoCap C18 Fast-Flow

## Physical Properties

Product Description	Bonded Phase	Meso-pore	End-capping	Max. Operating Pressure
MonoCap C18 Fast-flow	Octadecyl Groups	15 nm	Yes	22 MPa

Workable at a broad range of linear velocity from 0.5 to 5 mm/s without sacrificing efficiency and separation at high speed. The number of theoretical plates produced by MonoCap C18 Fast-Flow is nearly equivalent to a totally porous particle type capillary column packed with a 5 µm packing material. Columns are protected by either metal or PEEK hardware.

### Workable at High Flow Rates without Sacrificing Efficiency



**Conditions**

- Column : MonoCap C18 for Fast-flow (150 × 0.2 mm I.D.)
- Eluent : A) CH<sub>3</sub>CN B) H<sub>2</sub>O  
A/B=50/50 ,v/v
- Col.Temp. : Ambient
- Detection : UV 210 nm (MU701, Cell Volume 18 nL)
- Injection Vol. : 0.5 µL
- Sample : 1. Thiourea  
2. Acetophenone  
3. Benzene  
4. Toluene  
5. Naphthalene



## MonoCap C18 Nano-flow




### Physical Properties

Product Description	Bonded Phase	Meso-pore	End-capping	Max. Operating Pressure
MonoCap C18 Nano-flow	Octadecyl Groups	11 nm	Yes	22 MPa

MonoCap C18 Nano-flow produces higher number of theoretical plates compared to a totally porous particle type capillary column packed with a 3  $\mu\text{m}$  packing material. It can be operated at a wide range of flow rate with low back pressure and achieve very high sensitive results in Nano-LC-ESI/MS applications. Columns are protected by either metal or PEEK hardware.

## MonoCap C18 WideBore



### Physical Properties

Product Description	Bonded Phase	Meso-pore	End-capping	Max. Operating Pressure
MonoCap C18 WideBore	Octadecyl Groups	11 nm	Yes	22 MPa

The MonoCap C18 Fast-flow is also available in 0.5 mm I.D. size, which can be used at a wide range of flow rate from 6 to 100  $\mu\text{L}/\text{min}$  without sacrificing efficiency. The number of theoretical plates produced by MonoCap C18 WideBore is nearly equivalent to a totally porous particle type capillary column packed with a 5  $\mu\text{m}$  packing material. Columns are protected by a metal hardware.

## MonoCap C18 Trap Column

### Physical Properties

Product Description	Bonded Phase	Meso-pore	End-capping	Max. Operating Pressure
MonoCap C18 Trap Column	Octadecyl Groups	11 nm	Yes	20 MPa

MonoCap C18 Trap columns have a relatively big throughpore and workable at a high flow rate such as 10  $\mu\text{L}/\text{min}$ . This benefit makes MonoCap C18 Trap columns to be appropriate for on-line preconcentration or desalting of protein and peptide samples prior to HPLC separation with mass spectrometry detection. End-fittings are 1/16" (10-32 UNF). 1/32" end-fittings are also available upon request.



## MonoCap Amide



### Physical Properties

Product Description	Bonded Phase	Meso-pore	End-capping	Max. Operating Pressure
MonoCap Amide	Carbamoyl Groups	15 nm	None	22 MPa


Amide groups are chemically bonded to the monolithic silica and makes it suitable for the analysis of sugars via HILIC mode. As the back pressure is significantly low, a 500 mm length MonoCap Amide column deliver over 40,000 plates offering high efficiency. Generally, HILIC mode uses acetonitrile at a concentration between 65-95 % in an aqueous buffer such as ammonium acetate or ammonium formate, which have high solubility in organic solvents. Columns are protected by either metal or PEEK hardware.

## MonoCap SCX

### Physical Properties

Product Description	Bonded Phase	Meso-pore	End-capping	Max. Operating Pressure
MonoCap SCX	Benzenesulfonyl Groups	11 nm	None	20 MPa

MonoCap SCX is bonded with benzene sulfonic acid groups (strong cation exchange) and appropriate for 2D LC applications for the separation of biomolecules such as peptides and proteins.





## Ordering Information

### MonoCap C18 Fast-Flow

- For end-fittings information, please refer to page 30.
- All 50 mm length PEEK columns does not come with a hardware and will be supplied with 3 pcs of columns only.

Description	I.D. (mm)	Length (mm)	Hardware	Qty	Cat. No.
MonoCap C18 Fast-Flow	0.05	50		1/pk	5020-10102
		150	Metal	1/pk	5020-10101
		250		1/pk	5020-10100
		50		3/pk	5020-10002
		150	PEEK	1/pk	5020-10001
		250		1/pk	5020-10000
	0.075	50		1/pk	5020-10211
		150	Metal	1/pk	5020-10212
		250		1/pk	5020-10213
		50		3/pk	5020-10221
		150	PEEK	1/pk	5020-10222
		250		1/pk	5020-10223
	0.1	50		1/pk	5020-10112
		150	Metal	1/pk	5020-10111
		250		1/pk	5020-10110
		50		3/pk	5020-10012
		150	PEEK	1/pk	5020-10011
		250		1/pk	5020-10010
	0.2	50		1/pk	5020-10122
		150	Metal	1/pk	5020-10121
		250		1/pk	5020-10120
		50		3/pk	5020-10022
		150	PEEK	1/pk	5020-10021
		250		1/pk	5020-10020

**Ordering Information**

**MonoCap C18 Nano-Flow**

- For end-fittings information, please refer to page 30.
- All 50 mm length PEEK columns does not come with a hardware and will be supplied with 3 pcs of columns only.

Description	I.D. (mm)	Length (mm)	Hardware	Qty	Cat. No.
MonoCap C18 Nano-Flow	0.05	50	Metal	1/pk	5020-10143
		150		1/pk	5020-10141
		50	PEEK	3/pk	5020-10043
		150		1/pk	5020-10041
	0.075	50	Metal	1/pk	5020-10231
		150		1/pk	5020-10232
		50	PEEK	3/pk	5020-10241
		150		1/pk	5020-10242
	0.1	50	Metal	1/pk	5020-10153
		150		1/pk	5020-10151
		50	PEEK	3/pk	5020-10053
		150		1/pk	5020-10051
	0.2	50	Metal	1/pk	5020-10163
		150		1/pk	5020-10161
		50	PEEK	3/pk	5020-10063
		150		1/pk	5020-10061

**MonoCap C18 WideBore**

- For end-fittings information, please refer to page 30.

Description	I.D. (mm)	Length (mm)	Hardware	Qty	Cat. No.
MonoCap C18 WideBore	0.5	50	Metal only	1/pk	5020-10202
		150		1/pk	5020-10201
		250		1/pk	5020-10200

## Ordering Information

### MonoCap C18 Trap Column

- For end-fittings information, please refer to page 30.

Description	I.D. (mm)	Length (mm)	Hardware	Qty	Cat. No.
MonoCap C18 Trap Column	0.05	50	With Hardware	1/pk	5020-10026
		100		1/pk	5020-10038
		150		1/pk	NA
		50	Without Hardware	1/pk	5020-10027
		100		1/pk	5020-10039
		150		1/pk	NA
	0.075	50	With Hardware	1/pk	5020-10028
		100		1/pk	5020-10036
		150		1/pk	NA
		50	Without Hardware	1/pk	5020-10029
		100		1/pk	5020-10037
		150		1/pk	NA
	0.2	50	With Hardware	1/pk	5020-10033
		100		1/pk	NA
		150		1/pk	NA
50		Without Hardware	1/pk	5020-10034	
100			1/pk	NA	
150			1/pk	5020-10031	

### MonoCap Amide

- For end-fittings information, please refer to page 30.

Description	I.D. (mm)	Length (mm)	Hardware	Qty	Cat. No.
MonoCap Amide	0.075	150	Metal	1/pk	5020-10191
		250		1/pk	5020-10192
		500		1/pk	5020-10193
		150	PEEK	1/pk	5020-10091
		250		1/pk	5020-10092
		500		1/pk	5020-10093
	0.1	150	Metal	1/pk	5020-10181
		250		1/pk	5020-10182
		500		1/pk	5020-10183
		150	PEEK	1/pk	5020-10081
		250		1/pk	5020-10082
		500		1/pk	5020-10083
	0.2	150	Metal	1/pk	5020-10171
		250		1/pk	5020-10172
		500		1/pk	5020-10173
150		PEEK	1/pk	5020-10071	
250			1/pk	5020-10072	
500			1/pk	5020-10073	

**Ordering Information**

**MonoCap SCX**

- For end-fittings information, please refer to the following information.

Description	I.D. (mm)	Length (mm)	Hardware	Qty	Cat. No.
MonoCap SCX	0.2	50	Metal	1/pk	5020-10174
		150		1/pk	5020-10175
		250		1/pk	5020-10176
		500		1/pk	5020-10177
		50	PEEK	1/pk	5020-10074
		150		1/pk	5020-10075
		250		1/pk	5020-10076
		500		1/pk	5020-10077

**End-fittings of MonoCap Monolithic Capillary HPLC Columns**

Description	End-fittings Details
MonoCap C18 Fast-flow MonoCap C18 Nano-flow MonoCap C18 WideBore MonoCap Amide MonoCap SCX	1. Metal Hardware Type End-fittings are Valco 1/16" (10-32 UNF). Valco 1/32" (6-40 UNF) end-fittings can also be arranged upon request, indicate "1/32" when ordering.
	2. PEEK Hardware Type 1/16" male nut, ferrule and PTFE sleeve are included.

**Connection Kit for MonoCap C18 Trap Column**

Description	Cat. No.
MonoCap C18 Trap Column Connection Kit 1/16" (Union·Sleeve·Capillary Tubing 2 pcs each, Nut·Ferrule 4 pcs each)	5020-10044
MonoCap C18 Trap Column Connection Kit 1/32" (Union·Sleeve·Capillary Tubing 2 pcs each, Nut·Ferrule 4 pcs each)	5020-10045
MonoCap C18 Trap Column Assembly Parts 1/16" (Nut·Ferrule 4 pcs each)	5020-10046
MonoCap C18 Trap Column Assembly Parts 1/32" (Nut·Ferrule 4 pcs each)	5020-10047

