

Simultaneous determination of VMA, HVA, and creatinine

Neuroblastomas (NB) are catecholamine producing tumors and diagnosis is possible by analyzing for metabolic by-products of catecholamine in urine.

In this report vanillylmandelic acid (VMA), homovanillic acid (HVA) and creatinine were measured simultaneously by reversed phase HPLC and pre-processing was remarkably improved by HPLC autosampling with a micro-plate (Fig. 1).

Fig. 2 shows chromatograms of standard samples. A UV detector was used for creatinine and a coulometric detector (ECD) for VMA and HVA.

Fig. 3 shows the chromatogram from a urine test paper sample. Specialized data processing software is being prepared.

Conditions:

Column: Catechopak II
 Eluent: 50mMKH₂PO₄/ CH₃CN
 = 2745 / 255 EDTA 10mg/L
 1-Octanesulfonic acid Na 3mM
 pH 2.0(with H₃PO₄)
 Wave length(UV): 235nm
 Electrode potential: D1 0.2V, D2 0.25V, G 0.30V
 Flow rate: 1.2ml/min
 Column temperature: 45 degree celsius
 Injection volume: 10ul
 Sample: STD mixture
 VMA,VLA,HVA:250ppb
 Creatinine 25ppm




Fig. 1 Pre-preparation procedure




Fig. 2 Standard chromatogram




Fig. 3 Sample (6 months old urine test paper) chromatogram