

# How to Use ISOLUTE® PLD+ Protein and Phospholipid Depletion Plates

ISOLUTE® PLD+ plates provide extremely efficient removal of proteins and phospholipids from blood based samples using a simple 3 step procedure.

## Recommended Operating Procedure

Before you start, ensure collection plate is in position.

1. Dispense 300–400  $\mu\text{L}$  crash solvent into each well. No dripping will occur.
2. Dispense 100  $\mu\text{L}$  sample (e.g. serum, plasma, whole blood) into each well. To ensure full protein precipitation and avoid cloudy extracts, it is important to achieve efficient mixing on sample addition. This can be achieved by:
  - a. Vortex mix: 30s-1min, or
  - b. Aspirate/dispense: 3–4 cycles, or
  - c. Dispense sample vertically into well, with force, and wait 5 minutes
3. Apply vacuum or positive pressure until sample is fully eluted (5 mins).

### Recommended Crash Solvent

Acetonitrile or acetonitrile containing 1% (v/v) formic acid. Methanol or acidified methanol is also suitable.

See FAQ's ([www.biotage.com](http://www.biotage.com)) for further information on maximizing analyte recovery.

### Mixing

If insufficient mixing occurs on sample addition, proteins will not precipitate efficiently in the well. This can lead to turbidity (cloudiness) in extracts or well blockage during processing.

### Positive pressure processing:

Recommended pressure: 3 psi. When processing using Biotage® Extrahera™, set to 0.3–0.4 bar, for 6 minutes.

### Vacuum processing:

Recommended vacuum: -0.2 bar

For maximum phospholipid removal, do not exceed these vacuum/pressure conditions. For very viscous samples, a slight increase in pressure or vacuum may be required. The use of high vacuum/pressure conditions may lead to breakthrough of matrix components.




Figure 1. Typical ISOLUTE® PLD+ procedure.

## Ordering Information

Part Number	Description	Quantity
<b>918-0050-P01</b>	ISOLUTE® PLD+ Protein and Phospholipid Removal Plate	1
<b>Accessories</b>		
<b>121-5202</b>	Collection plate, 1 mL, square	50
<b>121-5203</b>	Collection plate, 2 mL, square	50
<b>121-5213</b>	Collection plate, 2 mL, round	50
<b>121-5204</b>	Piercable sealing cap	50
<b>Vacuum processing</b>		
<b>121-9600</b>	Biotage® VacMaster™-96 sample Processing manifold	1
<b>121-9601</b>	VacMaster VCU-1 Vacuum Control Unit	1
<b>121-9602</b>	VacMaster VCU-2 Vacuum Control and Generation Unit	1
<b>Positive pressure processing</b>		
<b>PPM-96</b>	Biotage® PRESSURE+ Positive Pressure Manifold, 96 position	1

### EUROPE

Main Office: +46 18 565900  
 Toll Free: +800 18 565710  
 Fax: +46 18 591922  
 Order Tel: +46 18 565710  
 Order Fax: +46 18 565705  
 order@biotage.com  
 Support Tel: +46 18 56 59 11  
 Support Fax: + 46 18 56 57 11  
 eu-1-pointsupport@biotage.com

### NORTH & LATIN AMERICA

Main Office: +1 704 654 4900  
 Toll Free: +1 800 446 4752  
 Fax: +1 704 654 4917  
 Order Tel: +1 704 654 4900  
 Order Fax: +1 434 296 8217  
 ordermailbox@biotage.com  
 Support Tel: +1 800 446 4752  
 Outside US: +1 704 654 4900  
 us-1-pointsupport@biotage.com

### JAPAN

Tel: +81 3 5627 3123  
 Fax: +81 3 5627 3121  
 jp\_order@biotage.com  
 jp-1-pointsupport@biotage.com

### CHINA

Tel: +86 21 2898 6655  
 Fax: +86 21 2898 6153  
 cn\_order@biotage.com  
 cn-1-pointsupport@biotage.com

To locate a distributor,  
 please visit our website at  
[www.biotage.com](http://www.biotage.com)

#### Part Number: U1314.V.1

© 2015 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. Product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.

For more information visit [www.biotage.com](http://www.biotage.com).