


# Extraction Cortisol From Human Saliva Using ISOLUTE® SLE+ Plates Prior to LC-MS/MS Analysis

This application note describes the supported liquid extraction clean-up of Cortisol from saliva prior to quantitative LC-MS/MS analysis.



**Figure 1.** Structure of Cortisol

## Introduction

Cortisol is a steroid hormone that when measured from saliva can be used as an indication of stress. This methodology has been designed to give an effective and efficient supported liquid extraction protocol for the clean-up and concentration of salivary cortisol levels.

Analyte recoveries achieved using this method ranged from 96-99% with RSDs below 3% (n=7) for all analytes.

ISOLUTE SLE+ Supported Liquid Extraction plates offer an efficient alternative to traditional liquid liquid extraction (LLE) for bioanalytical sample preparation, providing high analyte recoveries, no emulsion formation, and significantly reduced sample preparation time.

## Analytes

Cortisol, Cortisol D<sub>4</sub>.

## Sample Preparation Procedure

<b>Plate configuration:</b>	ISOLUTE SLE+ 200 µL Supported Liquid Extraction Plate, part number 820-0200-P01
<b>Sample pre-treatment:</b>	Dilute sample 1:1 (v/v) with water.
<b>Sample loading:</b>	Load the pre-treated sample (200 µL total volume) onto the plate and apply a pulse of vacuum (VacMaster 96 Sample Processing Manifold, 121-9600) or positive pressure (PRESSURE+ 96 Positive Pressure Manifold PPM-96) to initiate flow. Allow the sample to adsorb for 5 minutes.
<b>Analyte extraction:</b>	Apply MTBE (1 mL) and allow to flow under gravity for 5 minutes. Apply vacuum or positive pressure to pull through any remaining extraction solvent, collecting into collection plate.
<b>Post extraction:</b>	Evaporate the extract to dryness (40 °C) (SPE Dry 96 Sample Concentrator SD-9600-DHS-EU). Reconstitute in water:methanol (50:50, v/v) (100 µL).

## HPLC Conditions

**Instrument:** Waters Aquity UPLC

**Column:** Aquity UPLC BEH C18 1.7  $\mu\text{m}$  x 2.1 x 50 mm

**Mobile Phase:** A: 2 mM  $\text{NH}_4\text{OAc}$  0.1% formic acid (aq)  
B: 2 mM  $\text{NH}_4\text{OAc}$  0.1% formic acid in methanol

Gradient:	Time (mins)	%A	%B	Flow	Curve
	0.00	50	50	0.4	1
	0.8	50	50	0.4	6
	3.30	19	81	0.4	6
	3.50	50	50	0.4	1

**Injection:** 10  $\mu\text{L}$  (partial loop with needle overfill)

**Sample temperature:** 20  $^{\circ}\text{C}$

**Column temperature:** 40  $^{\circ}\text{C}$

## Mass Spectrometry Conditions

**Instrument:** Waters Quattro Premier XE triple quadrupole mass spectrometer equipped with an electrospray interface

**Source temperature:** 150  $^{\circ}\text{C}$


**Desolvation temperature:** 450  $^{\circ}\text{C}$

**Table 1.** MRM transitions in positive mode


Compound	MRM	Dwell (s)	Cone voltage	Collision energy
Cortisol quant	363.30 > 121.00	0.08	25	25
Cortisol qual 1	363.30 > 105.00	0.08	25	45
Cortisol qual 2	363.30 > 91.00	0.08	25	45
Cortisol D4 (ISTD)	367.30 > 121.00	0.08	25	25

## Results

This ISOLUTE SLE+ protocol demonstrates analyte recoveries ranges from 96-99% as shown in figure 3 (page 3) with RSDs below 3% for both Cortisol and Cortisol D4. Figure 2 shows a typical chromatogram at a concentration range of 25ng/mL, with a calibration curve demonstrating linearity from 25-2000 ng/mL.



**Figure 2.** Typical Chromatogram showing extracted Cortisol at 25 ng/mL and calibration curve showing linearity from 25-2000 ng/mL



**Figure 3.** Typical analyte % recoveries for Cortisol and Cortisol D4 (n=7) using the ISOLUTE SLE+ protocol

In accordance with known clinical detection requirements this assay was demonstrated at a level of 25 ng/mL, however based upon the achieved linearity and signal to noise ratio, it is estimated that a lower limit of quantitation of 0.5 ng/mL could be reached.

## Ordering information

Part Number	Description	Quantity
<b>820-0200-P01</b>	ISOLUTE SLE+ 200 µL Supported Liquid Extraction Plate	1
<b>PPM-96</b>	PRESSURE+96 Positive Pressure Manifold	1
<b>121-9600</b>	VacMaster 96 Sample Processing Manifold	1
<b>SD-9600-DHS-EU</b>	SPE Dry 96 Sample Concentrator	1

To search and download more of Biotage's extensive database of sample preparation applications please visit <http://www.biotage.com/applications> or scan the QR code with your smart phone to go direct.



**EUROPE**  
Main Office: +46 18 56 5900  
Fax: +46 18 59 1922  
Order Tel: +46 18 56 57 10  
Order Fax: +46 18 56 57 05  
order@biotage.com  
EU-1pointsupport@biotage.com

**NORTH AMERICA**  
Main Office: +1 704 654 4900  
Toll Free: +1 800 446 4752  
Fax: +1 704 654 4917  
Order Tel: +1 704 654 4900  
Order Fax: +1 434 296 8217  
ordermailbox@biotage.com  
US-1pointsupport@biotage.com

**JAPAN**  
Tel: +81 3 5627 3123  
Fax: +81 3 5627 3121  
jp\_order@biotage.com  
JP-1pointsupport@biotage.com

**CHINA**  
Tel: +86 21 2898 6655  
Fax: +86 21 2898 6153  
cn\_order@biotage.com  
CN-1pointsupport@biotage.com  
www.biotage.com

For regional distributors, products and offers, please visit [www.biotage.com](http://www.biotage.com)