


# Extraction of Amphetamines From Urine Using ISOLUTE® SLE+ 96-Well Plates

## Introduction

This application note describes the extraction of amphetamines from urine using ISOLUTE SLE+ 96 fixed-well plates with LC-MS/MS analysis.



**Figure 1.** Structure of Amphetamine

This method describes the use of ISOLUTE SLE+ in 96 well plate format for the extraction of a range of amphetamines. This procedure is effective with sample volumes of 100 µL for analyte recoveries > 90%. For larger volume (500 µL) samples, see AN746.

ISOLUTE SLE+ Supported Liquid Extraction plates and columns offer an efficient alternative to traditional liquid-liquid extraction (LLE) for bioanalytical sample preparation, providing high analyte recoveries, no emulsion formation, and significantly reduced sample preparation time.

## Analytes

Ephedrine, amphetamine, methamphetamine, MDMA, MDA, MDEA.

## ISOLUTE SLE+ procedure

ISOLUTE SLE+ 200 Supported Liquid Extraction Plate, part number 820-0200-P01.

- Sample Pre-treatment:** Dilute urine (100 µL) 1:1 with 0.5 M ammonium hydroxide (100 µL).
- Sample Load:** Load pre-treated sample (200 µL) to plate followed by a pulse of vacuum to initiate flow and leave for five minutes. Add 0.05 M HCl/Methanol (100µL) to each collection plate well. \*
- Analyte Elution:** Elute with dichloromethane (1 mL). Leave to flow under gravity for 5 minutes then apply a second aliquot of dichloromethane (1 mL) followed by a short pulse of vacuum.
- Post extraction:** Evaporate to dryness and reconstitute in 200 µL 0.1% formic acid (aq) in water/methanol (90/10, v/v).
- Additional information:** \*During evaporation, amphetamines in the free base form suffer major losses. As a result it was necessary to convert these to a stable HCl form at low pH using 0.05M HCl/Methanol (100uL) in the collection plate. All samples were processed and dried down using the Vacmaster and SPE Dry respectively.

## UPLC Conditions

**Instrument:** Waters Acquity UPLC interfaced to a Quattro Premier XE triple quadrupole MS using electrospray ionization.

**Column:** Acquity BEH C18 100 x 2.1mm x 1.7u.

**Mobile phase:** Isocratic, 0.1% formic acid aq/0.1% formic acid in methanol (80/20, v/v).

**Flow rate:** 0.43 mL/min.

**Temperature:** 40 °C.

## Mass Spectrometry Conditions

**Source temp:** 150 °C.

**Desolvation temp:** 450 °C.


**Collision cell pressure:** 3.58 e<sup>-3</sup> mbar.

**Table 1.** MRM transitions for a range of amphetamines.


Scan function	Compound	MRM transition	Cone voltage	Collision energy
1	Ephedrine	166.1 > 133.0	20	19
2	Amphetamine	136.0 > 118.9	16	9
	Methamphetamine	150.0 > 90.9	22	17
	MDA	180.1 > 105.0	16	23
	MDMA	194.1 > 163.0	20	13
3	MDEA	208.2 > 163.0	22	13

## Results

Figure 2 shows the total ion chromatogram for all amphetamine analytes detected using this methodology whilst recoveries of all analytes were > 90% at an LLOQ of 500 pg/mL using the standard SLE+ 200 plate method as shown in figure 3.



**Figure 2.** Total ion chromatogram for all amphetamine analytes at the LLOQ level of 500 pg/mL



**Figure 3.** Analyte % recoveries of amphetamines using SLE+ 200 methodology at 500 pg/mL (n=7)

### Ordering information

Part number	Description	Quantity
820-0200-P01	ISOLUTE SLE+ 200 µL supported liquid extraction plate	1
121-9600	VacMaster-96 Sample processing manifold complete (without Vacuum control)	1
SD2-9600-DHS-UK	SPE Dry 96 Dual, 240 V UK	1

#### NORTH AMERICA

Main Office: +1 704 654 4900  
Toll Free: +1 800 446 4752  
Fax: +1 704 654 4917  
Order Tel: +1 704 654 4900  
press (4) at the auto attendant  
Order Fax: +1 434 296 8217  
ordermailbox@biotage.com  
1-pointsupport@biotage.com

#### EUROPE

Main Office: +46 18 56 5900  
Fax: +46 18 59 1922  
Order Tel: +46 18 56 57 10  
Order Fax: +46 18 56 57 05  
order@eu.biotage.com

#### Japan

Tel: +81 3 5627 3123  
Fax: +81 3 5627 3121  
jp\_order@biotage.com

#### Distributors

To locate a distributor please visit our  
Web site at [www.biotage.com](http://www.biotage.com)