

Multi-wavelength UV/VIS Detector

UV-2077



The Model UV-2077 is a multi-wavelength UV/VIS detector based upon a photo-diode array detection system. The state-of-the-art optics and electronics are housed in the same compact case as other LC-2000Plus series detectors including the UV-2070. The UV-2077 has a unique design feature that separates it from both conventional detectors which employ scanning/readout shift registers or incorporates a 36-element discrete diode-array, enabling truly simultaneous 36-channel detection in real-time with maximum sensitivity and selectivity. This configuration allows simultaneous monitoring of chromatograms at 4 different wavelengths or 3 wavelengths and a ratio of 2 others. The advanced computer technology allows the wavelengths to be selected in 1nm increments. In addition the detector is able to acquire spectra truly on-the-fly without interrupting the chromatogram. The UV-2077 allows the analyst to obtain the high quality chromatographic and spectral data of an expensive PDA detector in an affordable, compact and easy to use package.

The UV-2077 can easily be interfaced with the JASCO HPLC system using the new 32 bit HSS-2000 HPLC system software for Windows2000. The HSS-2000 package offers full automation of the JASCO HPLC system including the UV-2077.

- Significantly lowered baseline noise and drift due to unique electronics and optics
- Uninterrupted chromatographic peaks and on-flow spectral collection utilizing single beam photometry
- Compact 150 (W) x 150 (H) x 470 (D)mm size matching other LC-2000Plus Series modules

SPECIFICATIONS

Optical system

Photometric system: Flat field type polychromator with concave grating
Wavelength range: 200 ~ 600 nm, settable in 1 nm increment
Wavelength accuracy: ± 7 nm
Light source: Deuterium lamp

Photometric system

Photometric method: Single beam photometry
Detector: 36 elements photodiode
Noise level: $\pm 1.0 \times 10^{-5}$ ABU at 250 nm with empty cell (based on ASTM method)

Drift: $\pm 5 \times 10^{-4}$ ABU/h
Response: 3 steps (FAST, STD and SLOW)

Flow through cell

Mounting: Front loading cassette type with bulkhead union
Flow cell: Tapered cell

Signal input/output

Analog outputs: 4 channel independent outputs; from no. 1 to no. 3 channels for absorbance output, no. 4 channel; absorbance or ratio output, selectable

Output to recorder: 10 mV/full scale
Output to integrator: 1V/AU, 0.2V/0.2AU, 1V/2AU
Others: Maker output/input, Autozero input, LC-NET output/input, Program reset/run, Leak out

Self-diagnostics

For starting up: START UP, RAM, ROM, ± 5 V, ± 12 V, ADC CHK RUNNING, LEAK IN CELL, D2 LAMP OFF, TIME OFF

Program function

FTime program: Wavelength, range, auto zero, baseline correction, response, spectrum measurement (0.1 min increment, up to 999.9 min, 64 steps)
Spectral memory: up to 32 spectra can be stored
Auto lamp off: Settable in 0 ~ 99 hour, 0.1 hour increment

Others

Dimensions: 150 (W) x 470 (D) x 150 (H) mm
Weight: Approx. 9 kg
Power Input voltage: AC 100 ~ 240 V $\pm 10\%$, 50/60Hz
Power requirements: 100VA

ORDERING INFORMATION

6826-J002A	Model UV-2077	Multi-wavelength UV/VIS Detector
6824-J002A	Model UV-2070	UV/VIS Detector (~ 900 nm)
6825-J002A	Model UV-2075	UV/VIS Detector (~ 600 nm)
6832-J002A	Model MD-2015	Diode Array Detector (~ 900 nm)
6831-J002A	Model MD-2010	Diode Array Detector (~ 600 nm)
6829-J003A	Model FP-2020	Fluorescence Detector
6833-J002A	Model RI-2031	Refractive Index detector
6838-J002A	Model CD-2095	Circular Dichroism Detector
6837-J002A	Model OR-2090	Chiral Detector
6836-J002A	Model CL-2027	Chemiluminescence Detector
6818-J002A	Model PU-2080	HPLC Pump
6823-J003A	Model PU-2089	Quaternary HPLC Pump
6821-J002A	Model PU-2086	Semi-preparative HPLC Pump
6822-J002A	Model PU-2087	Preparative HPLC Pump
6819-J002A	Model PU-2080i	Bio-inert HPLC Pump
6823-J004A	Model PU-2089i	Bio-inert Quaternary HPLC Pump
6820-J002A	Model PU-2085	Semi-micro HPLC Pump
6828-J008A	Model AS-2059	Autosampler
6827-J005A	Model AS-2055	Autosampler
6827-J006A	Model AS-2057	Autosampler
6834-J002A	Model CO-2060	Column Oven
6835-J002A	Model CO-2065	Column Oven

STANDARD ACCESSORIES

Fuse 2 pcs
AC power cable 1 pc



● Specifications are subject to change without notice.



JASCO INTERNATIONAL CO., LTD.

4-21, Sennin-cho 2-chome, Hachioji, Tokyo 193-0835, Japan Tel: +81-426-66-1322 Fax: +81-426-65-6512 Internet: <http://www.jascoint.co.jp/english/index.html>

JASCO INCORPORATED, 8649 Commerce Drive, Easton, Maryland 21601-9903, U.S.A Tel: +1-800-333-5272 Tel: +1-410-822-1220 Fax: +1-410-822-7526 Internet: <http://www.jascoinc.com>

JASCO EUROPE s.r.l., Via Confalonieri 25, 22060 Cremella (Co), Italy Tel: +39-039-956439 Fax: +39-039-958642 Internet: <http://www.jasco-europe.com>

JASCO Labor-und Datentechnik GmbH Deutschland, Robert-Bosch-Strasse 11, 64823 Gross-Umstadt, Germany Tel: +49-6078-74949 Fax: +49-6078-74006 Internet: <http://www.jasco.de>

JASCO (UK) Ltd., 18 Oak Industrial Park, Chelmsford Road, Great Dunmow, Essex, CM6 1XN, U.K. Tel: +44-1371-876988 Fax: +44-1371-875597 Internet: <http://www.jasco.co.uk>

JASCO Eastern Europe Information Office, H-1113, Karolina UT 29-31, Budapest, Hungary Tel: +36-1-4669549 Fax: +361-2093538

JASCO France s.a.r.l., 60 Bd du Maréchal Juin, 44100-Nantes, France Tel: +33-2-40430043 Fax: +33-2-40432437 Internet: <http://www.jascofrance.fr>

JASCO Benelux B.V., Industrieweg 68-3606 AS Maarssen, Netherlands Tel: +31-346-556699 Fax: +31-346-554619 Internet: <http://www.jasco.nl>

JASCO Analitica Spain S.L., Emiliano Barral 13-B, 28043 Madrid, Spain Tel: +34-91-361-2440 Fax: +34-91-356-7058 Internet: <http://www.jasco-spain.com>

JASCO Scandinavia AB, Ågatan 32, 431 35 Mölndal, Sweden Tel: +46-31-27 03 60 Fax: +46-31-746 85 90 Internet: <http://www.jascoscandinavia.se>

Serving the global marketplace
with Analytical Instrumentation

JASCO Corporation
2007-16, Mikasa-cho, Hachioji, Tokyo 193-0837, Japan
<http://www.jasco.co.jp>