3P sim



A modeling and evaluation tool for dynamic sorption data

TOOLS FOR DISCOVERY



Characterization of

particles • powders • pores

info@3P-instruments.com www.3P-instruments.com www.dynamicsorption.com

PORE DETERMINATION



for dynamic sorption data

Introduction

Modeling of dynamic sorption processes are very useful for the understanding of sorption characteristics on fixed bed adsorbers. With the help of such tools, the experimental time can be drastically reduced and parametric studies can be easily performed. The simulation and evaluation tool **3P sim** is provided with the new dynamic sorption analyzer series "mixSorb" for the investigation of breakthrough characteristics. This software allows the calculation of equilibrium data of mixtures from pure component isotherms as well as the evaluation of breakthrough curves, based on mass- and energy balances. The available solver is easy to handle and can be operated by untrained operators. As an example, the calculation of CH_4 / CO_2 breakthrough curves on activated carbon is shown.



Gas Purification and Separation



Figure:

Different calculations with mixture models using a 3D-plot **black:** total loading green: partial loading component A **blue:** partial loading component B **blue labelling:** calculation with constant gas phase composition red labelling: calculation with constant total pressure



Isotherm Fitting

With **3P sim** measured isotherm data can be fitted with the following mathematical isotherm models:

- HENRY
- LANGMUIR
- TOTH
- SIPS
- FREUNDLICH
- DUALSITE LANGMUIR
- DUALSITE LANGMUIR SIPS



The program allows the calculation of total and partial loadings at given pressures or compositions and supports the following theories:

- IAST* with LANGMUIR
- IAST with TOTH
- IAST with DUALSITE LANGMUIR
- IAST with DUALS. LANGMUIR SIPS
- Multicomponent LANGMUIR
- Multicomponent SIPS

* Ideal Adsorption Solution Theory

Dynamic Simulation

- 3P sim provides solutions for massand energy balances which allow simulations without user precognition or programming skills.
- Technically relevant transport parameters (e.g., LDF* constants) are accessible
- Simulation of breakthrough curves and temperature profiles

* Linear Driving Force





Mathematical description of the Isotherms for further calculation. With a set of measured isotherms, the temperature dependence can be calculated and a theoretical interpolated isotherm at a given temperature can be predicted.





Mixture equilibria can be implemented in Dynamic Simulations.





Complete data set for **parametric** studies.

3P sim

Conclusion

By measuring breakthrough curves of mixtures in a carrier gas, separation effects can be observed. If the isotherms are well known, a fitting of such breakthrough curves with a dynamic model is possible to get information for further parametric studies. With the dynamic method, it is possible to investigate adsorbents regarding their separation performance under relevant conditions. Such valuable performance characteristics can only be derived from dynamic experiments. **3P sim** allows parametric studies to reduce experimental effort and improves the understanding of separation processes.



Dynamic Sorption Analyzer

Industrial adsorbents such as active carbons, zeolites and silica gels are widely used in adsorptive separation processes on a multi-ton scale. The after-treatment of exhaust gasses, the removal of carbon dioxide in bio gas plants, purification and fractionation of natural gas, air separation, respiratory protection and separation of isomers are just a few examples where adsorptive separation is employed as the most efficient and economic separation technique.

Features and Benefits

Large Flow Range

Investigation of mixture equilibrium and kinetics, downscaling of technical processes.

 Built-in Thermal Conductivity Detector (TCD)

Automated time-resolved measurement of outlet gas composition.

- Built-in Sample Preparation up to 400°C Investigation of hydrophilic materials (i.e., zeolites and silica gels).
- Automated Pressure Regulation Completely programmable for operator-free performing of pressure steps.
- Automated Gas and Vapor Mixing No need of multiple pre-mixed gas tanks, increasing of flexibility of measurement conditions.
- Optional Gas Analysis via Interfaced MS Investigation of ternary and more complex systems.

Specifications

World-class performance is ensured by designing to very stringent instrument specifications thus ensuring the highest quality data.

Adsorber:	1
No. of Mass Flow Controllers*:	Up to 4
Max. Pressure:	Up to 68 bar
Sample Preparation:	In-situ
Temperature Range:	-20 – 400 °C
Height:	960 mm (38 in)
Width:	860 mm (34 in)
Depth:	640 mm (25 in)
Weight:	120 kg (265 lbs)

^{*} Mass flow controllers are available in different ranges. For optimal instrument configuration with respect to customer's needs.



YOUR PARTNER IN PARTICLE CHARACTERIZATION

3P Instruments offers profound expertise for the characterization of particles, powders and pores in Europe.



www.dynamicsorption.com

collects and summarizes all the features, advantages, examples and knowledge of dynamic sorption (or flow sorption) methods. Feel free to browse through the pages about breakthrough curves, mixture adsorption and the commercial breakthrough analyzers of the mixSorb series, and our specialist literature for download under "Resources".



Characterization of particles • powders • pores

info@3P-instruments.com www.3P-instruments.com Tel. +49 8134 9324 0 www.3P-instruments.com