

Higher-Order Structure Analysis of High Concentration Monoclonal Antibody by Circular Dichroism (CD) and Infrared (IR) Spectroscopy

Ai Yamane¹, Taiji Oyama¹, Kohei Tamura¹, Miyuki Kanno¹, Shingo Norimoto¹, Forrest Kohl², Satoko Suzuki¹, and Kenichi Akao¹

¹JASCO Corporation, Hachioji, Tokyo, 192-8537, Japan, ²JASCO Inc., Easton, MD, USA

fkohl@jascoinc.com

INTRODUCTION Therapeutic antibodies have been dramatically expanding their market over the past decade and become one of the major therapeutic proteins. Among different physical properties of therapeutic antibodies, higher-order structure (HOS) is an integral part of their characterization¹. One of the most well-known methods for HOS characterization is circular dichroism (CD) spectroscopy. It provides both secondary and tertiary structure information of relatively low-concentrated protein solution (< ≈10 mg/mL). Infrared (IR) spectroscopy is also a widely used method to analyze the secondary structure of proteins. It can measure not only high-concentrated protein solutions (100 mg/mL or higher), but also suspensions and solids. As different techniques for HOS analysis exist, orthogonal method using more than one technique has been suggested for full characterization of HOS^{1,2}. In this study, we estimated the secondary structure contents of rabbit serum IgG at concentrations of 5.8 to 23.1 mg/mL from its far-UV CD and IR spectrum. For CD secondary structure analysis (SSE), we used the J-1500 CD spectrometer which can measure CD and absorbance simultaneously, and the BeStSel³ algorithm which shows high accuracy of β -sheet estimation. For IR SSE, we used our latest FTIR spectrometer, FT/IR-4X, and our lately improved IR SSE software, IR SSE-4000. Furthermore, we also succeeded in acquiring near-UV CD spectrum of 5.8 to 156.2 mg/mL of IgG by selecting appropriate pathlengths for each concentration.

EXPERIMENTAL

Sample

- IgG, from rabbit serum (Sigma Aldrich):
5.8 to 156.2 mg/mL
- Buffer: 20 mM citric acid buffer
(150 mM NaCl, pH 6.0)




Fig 1. Image from the RCSB PDB (<http://www.rcsb.org/>) of PDB ID 1HZH⁴.

Measurement and Analysis

CD Spectrum

J-1500 CD Spectrometer



Short pathlength cuvette



- CD/absorbance simultaneous measurement with high accuracy
- Proteins at a wide range of concentrations measured with the use of short pathlength cuvette

- Well-known CD SSE algorithm
- Highly accurate estimation of β -sheet content

IR Spectrum

FT/IR-4X FTIR Spectrometer

NEW



ATR PRO 4X Single-Reflection ATR

- High S/N and resolution with small body
- Proteins easily measured by putting one-drop of sample on the ATR PRO 4X

IR SSE-4X Secondary Structure Estimation Program

NEW

- SSE accuracy improved by the addition of protein spectra to the reference model

1. Secondary Structure Estimation by Far-UV CD and IR Spectrum

Sample concentrations: 5.8, 11.4, 17.5, and 23.1 mg/mL

Table 1. Measurement conditions for far-UV CD and IR spectrum.

Measurement	Pathlength/mm	Sample Volume/ μ L
Far-UV CD	0.01	1.8
IR	-	10

2. Near-UV CD spectrum measurement

Sample concentrations: 5.8, 11.4, 17.5, 23.1, 43.3, 61.7, 76.8, 114.9, and 156.2 mg/mL

Table 2. Measurement conditions for near-UV CD spectrum.

Measurement	Concentration/mg \cdot mL ⁻¹	Pathlength/mm	Sample Volume/ μ L
Near-UV CD	5.8 and 11.4	1	280
	17.5 and 23.1	0.5	220
	43.3 or higher	0.1	22

RESULTS

1. Secondary Structure Estimation by Far-UV CD and IR Spectrum




Fig 2. (a) CD spectra, (b) absorbance spectra, (c) MRE spectra at far-UV region, and (d) IR spectra of rabbit serum IgG.

- ✓ The concentration of each sample was calculated from the absorbance at 205 nm^{5,6} (b), then the CD spectra (a) were converted into MRE spectra (c) using the calculated concentration. The MRE spectrum did not change its shape despite the change in the concentrations.
- ✓ The IR spectrum (d) did not seem to have a significant change in its shape regardless of the concentrations.

Secondary structure Estimation (SSE)




Fig 3. Secondary structure contents estimated from (a) MRE spectra and (b) IR spectra. Secondary structure contents were estimated from the MRE spectrum using the BeStSel program, and from the Amide I band (1700 to 1600 cm⁻¹) on IR spectra using the IR SSE-4000.

- ✓ Secondary structure content did not show significant differences among these concentrations.
- ✓ The results of both CD and IR SSE showed very close secondary structure contents.

REFERENCES

- 1) US Food and Drug Administration. Regulatory Consideration for the Characterization of HOS in Biotechnology Products. 2016
- 2) WHO. Guidelines on the quality, safety, and efficacy of biotechnological protein products prepared by recombinant DNA technology. 2013
- 3) A. Micsonai, F. Wien, L. Kernya, Y.H. Lee, Y. Goto, M. Refregiers, J. Kardos, Accurate secondary structure prediction and fold recognition for circular dichroism spectroscopy, PNAS.112 (2015) 3095-3103.
- 4) E.O. Saphire, P.W. Parren, R. Pantophlet, M.B. Zwick, G.M. Morris, P.M. Rudd, R.A. Dwek, R.L. Stanfield, D.R. Burton, I.A. Wilson, Science, 2001, 293, 1155-1159.
- 5) Nicholas J. Anthis, G. Marius Clore, Protein Science, 2013, 22, 851-858.
- 6) A. Micsonai, E. Buluyaki, J. Kardos, Methods in Molecular Biology, 2020, 2199, 175-189.

2. Near-UV CD Spectrum Measurement




Fig 4. Normalized (a) CD spectra, and (d) absorbance spectra on near-UV region of rabbit serum IgG. Since each sample was measured with different pathlength, the measured CD and absorbance spectrum was normalized as if all the samples with different concentrations were measured with the same pathlength (0.1 mm).

- ✓ CD signals at 255 nm to 305 nm under both concentrations of (b) 5.8 mg/mL and (c) 156.2 mg/mL were well-detected, showing that IgG at a wide range of concentrations was successfully measured with good sensitivity at near-UV region.

CONCLUSIONS

- The agreement on CD and IR SSE results suggests their use on orthogonal analysis, or one of them can be selected properly. For example, the amide I band on IR spectrum can be overlapped with some type of buffer (e.g. citrate buffer), therefore CD is a better option in this case. On the other hand, while IR can handle protein with concentrations of 100 mg/mL or higher, those for far-UV CD measurement is limited up to several tens of mg/mL.
- Selecting an appropriate pathlength depending on the protein concentration enables near-UV CD measurement at a wide range of protein concentrations.

ACKNOWLEDGEMENT

We are grateful to Prof. Kardos for useful discussions about secondary structure estimation of proteins and the use of the BeStSel program. (BeStSel URL: <http://bestsel.elte.hu/index.php>)