

TURBISCAN[®] TRILAB – TECHNICAL SHEET

1) Technical specifications and recommendations

<i>Technology</i>	Static Multiple Light Scattering (SMLS)
<i>Wavelength (nm)</i>	LED 880nm
<i>Detectors</i>	Backscattering and Transmission (for further information, see technical application note “Description of SMLS method for the evaluation of mean particle size in concentrated dispersions. “)
<i>Displacement interval max. resolution</i>	5µm
<i>Maximum displacement velocity</i>	15mm/s
<i>Scan step resolution</i>	20µm
<i>Sample preparation</i>	Native sample directly in the measurement cell, no dilution, no destructive method
<i>Measurement cells volumes (proposed by Formulation)</i>	4 and 30mL Note 1 : 4mL cells requires an adaptor provided by Formulation Note 2 : Formulation can provide custom adaptors on demand
<i>Minimum sample volume</i>	1.5 mL
<i>Maximum sample volume</i>	30 mL
<i>Particle size range</i>	10nm – 1mm (depending on physical parameters)
<i>Concentration range</i>	10 ⁻⁴ % - 95% v/v (depending on physical parameters)
<i>Available calculation</i>	TSI, TSI bottom, TSI middle, TSI Top, TSI manual, videos, mean value, phase thickness (mm), migration rates, mean particle size d _{SMLS} , Hydrodynamic particle size d _H
<i>Sample typology</i>	<ul style="list-style-type: none"> ✓ Non absorbing samples (diffusive samples) : stability + d_{SMLS} ✓ Absorbing samples (black samples) : stability + d_H
<i>Particle size validity conditions</i>	<ul style="list-style-type: none"> ✓ For diameter based on optical parameters (d_{SMLS}), the determination is only valid for diffusive samples. ✓ d_{SMLS} available only with 30mL cells and for Turbiscan LAB expert (Turbiscan LAB Standard and Thermo are not calibrated for d_{SMLS} measurements)
<i>ISO compliant</i>	TR 13097/TR 18811/TR 13014
<i>Dimension & weight</i>	35 x 48 x 58 cm – 35kg
<i>Space required</i>	Keep a distance of 30 cm minimum from each side
<i>Recommended computer configuration</i>	<ul style="list-style-type: none"> ✓ 2.5 GHz processor ✓ 4 GB RAM ✓ SSD or 7200rpm HDD ✓ Minimum 1 USB Port ✓ Microsoft WINDOWS 10 ✓ Screen resolution 1920x1080
<i>Barcode reader</i>	Yes

<i>Primary calibration</i>	All instruments are factory calibrated (backscattering and transmission signals) based on specific and traceable standards. The instrument is provided with calibration check-up standards for onsite calibration control (see above). In case of signal intensity deviation, only trained Formulation's representatives can recalibrate the instrument
<i>Calibration standards for instrument check-up</i>	<p>Teflon (PTFE) for Backscattering and Silicon Oil for Transmission (provided in a black box and included with the instrument)</p> <p>Validity: Calibration standards set is valid for 10 years if kept in the black plastic case at standard room temperature and humidity (serial number two first digits indicate production year)</p> <p>Traceability: Calibration standards set is tracked by Formulation with supplier batch number identification</p> <p>Data sheet: available under request.</p>
<i>Calibration check-up frequency recommendation</i>	Instrument calibration check every 3 months
« what corresponds to the 00 :00 :00 scan time »	“it corresponds to the starting time of the very first scan”
<i>Maximum scan frequency</i>	25 seconds for one scan measurement in position 1
<i>Maximum Scan number</i>	600 scans
<i>Reproducibility (corresponds to manual reproducibility with 10%-1µm latex sample for 10 scans with sample in and out between 2 scans)</i>	0.05%
<i>Repeatability (corresponds to automatic reproducibility with 10% 1µm latex sample for 10 scans with sample remaining inside the instrument between 2 scans)</i>	0.05%
<i>Operating temperature</i>	We recommend to keep a constant temperature during scan acquisition. Temperature variation affects refractive indices and signals intensity
<i>Temperature range</i>	[20°C-60°C]
<i>Temperature adjustment step</i>	0.1°C
<i>Temperature accuracy</i>	Thermoprobe Pt 100 ceramic CEI 60751 class A 3 strands
<i>Thermal security</i>	Thermostat switch with switch off at 110°C Switch action : Trip free manual reset: UL M2 class rating
<i>Voltage</i>	Acceptable voltages are between 100V and 240V AC
<i>Power consumption</i>	900W
<i>Electromagnetic compatibility and safety</i>	See Declaration of conformity

<i>Radiation</i>	No laser, LED power<1mW, class 1 (non directly visible)
<i>Spare parts availability</i>	Spare parts availability is guaranteed for 10 years.
<i>Preventive maintenance</i>	Formulation recommends preventive maintenance every year (replacement of wearing parts, cleaning and calibration)

2) Consumables

Items	Consumables	Composition	Thermal stability
<i>Cells 20mL</i>	Glass cell	Flat bottom borosilicated glass superior quality 1 st hydrolytic class (type1) (external diameter 27.5x72.5mm height with wall thickness 1.2mm) Cells are disposable. Rinse-off may lead to measurement quality issues. Cells must be changed after 3 washing cycles maximum to ensure glass quality	700°C
	Cap	Polypropylene	Up to 130°C
	Seal	Butyl/PTFE	-40 to 120°C
<i>Adapter</i>	For 4mL cells	Anodized aluminum	600°C
	For crimp cell	Anodized aluminum	600°C
<i>Sampling rack</i>		Thermoplastic (POM)	-40 to 115°C
<i>Cells 4mL</i>	Glass cell	borosilicated glass superior quality 1 st hydrolytic class (type1) (external diameter 15mmx45 mm height)	700°C
	Cap	Phenolic cap	Very good resistance
	Seal	PTFE	-200°C to 260°C
<i>Viscous product & foam cells</i>	Glass cell	Borosilicated glass superior quality 1 st hydrolytic class (type1)	700°C
	Bottom stopper	PTFE	-200°C to 260°C
	O ring	FKM Fluorocarbon rubber	-20 to 200°C
	Cap	Polypropylene	Up to 130°C
	Seal	Butyl/PTFE	-40 to 120°C
<i>Crimp cells</i>	Glass cell	Flat bottom borosilicated glass cells superior quality 1 st hydrolytic class (type1) (22.5x75.5mm)	700°C Pressure: Up to 10 bars with adapted caps
	Stoppers or crimp caps	Not provided	NA

Note 1: Elements autoclavable separately, not guaranteed on closed cell

Note 2: for chemical compatibility, please consult compatibility data table

3) References

- **Consumable's listing**
- **IQ-OQ**
- **Declaration of conformity in the instrument user guide**
- **Backscattering and Transmission standards data sheets**