

LIGHTIGO FireFly

Workhorse
for fast elemental analysis

What is FireFly?

FireFly is an analytical instrument for rapid **multi-elemental analysis** and **chemical imaging**.

FireFly utilises the Laser-Induced Breakdown Spectroscopy (LIBS) - one of today's most promising technique of elemental analysis.

With FireFly, you can easily apply the key benefits of LIBS technology in your elemental analysis application.

Applications

LIBS is fast and flexible, sensitive to the most of the chemical elements with detection limits of 1-100 ppm. Potential applications are:

- Metallic, geology, plastic or biology materials
- Elemental imaging (mapping)
- Depth profiling of multiplayer materials
- Toxic / heavy metals detection
- Markers and nanoparticles detection



Equipment

- Motorised 3-axis stage, resolution: <1 μm
- Microscopy view with objective turret
- Motorised laser spot size changing
- Motorised wavelength switching
- Gas purge and extraction
- Up to 100 Hz repetition rate

Geology




Biology



Plastics



Metals



FireFly specification

Table contains a list of generally available system parameters. Individual values may differ depended on the particular configuration.

Motorised Stage	Travel	100 × 100 × 70 (X × Y × Z) mm.
	Parameters	<1 µm micro step / 5 µm full step, speed up to 20 mm/s.
Sample handling	Sample size	Max 100 × 100 × 70 mm (irregular shape allowed).
	Sample holders	For general use: Universal clamp holder. For standard pellets: 12 × 12 mm 2 × 30 mm 1 × 50 mm.
Sample View	Sample view	Plan Apo objectives on motorized objective turret
	Stage view	Varifocal objective, Global shutter +Trigger-In.
Laser Focusing	Objective	Macro lens, long depth of focus, f = 75 mm, AR for 266/532/1064 nm. Aspherical micro-focus lens, f = 20 mm, AR for 266/532/1064 nm.
	Spot size	Motorised in a range of ~ 20–150 µm.
Laser	Type	Diode-Pumped Solid State, 5 ns, up to 100 Hz
	Wavelength	1064 nm (60 mJ) 532 nm (30 mJ) 266 nm (12 mJ)
	Features	Motorised wavelength switching and energy regulation.
Detection	Spectrometers	Czerny-Turner or Echelle; 180-1100 nm, resolution < 0.05 nm.
	Detectors	CCD/CMOS or intensified detectors, up to 100 Hz full frame, <10 ns gate width.
Synchronisation	Output	TTL output for synchronizing external devices, precision < 10 ns.
Gas Modules	Gas Purge	Continuous or triggered inert gas purging (gas input required).
	Gas Extraction	Active gas/dust extraction system with changeable filters.
Accessories	Calibration lamp	Hg-Ar line lamp (independent unit).
Control	Hardware	PC, LCD, keyboard, mouse.
	Software	Control software + Lightigo ImageLab SW Interactive spectra / chemical map browser, user profiles, multispot / lin- escan / mapping modes, measurement guides, elements database: NIST + LIBS, chemometrics, etc.
Technical Data	Mains Supply	110-240V 50/60Hz 10A.
	Protection	Interlock system on doors, laser filters on windows.
	Working environment	Indoor, 10-30 °C.
Dimensions and Weight	W × D × H	720 × 940 × 1500 mm.
	Weight	300 kg.