



CAPSULE SINKERS

Distek Capsule Sinkers come in a wide variety of shapes and sizes to fit all your Dissolution Applications. The chart below highlights our complete offering including; Image, Part Number, Description, and Dimensions.

| | PART NUMBER | DESCRIPTION | LENGTH | DIAMETER |
|---|-------------|---|---------|----------|
|  | 0500-0467 | Sinker, Basket with Clip on Lid, 316 SS | 26.9 mm | 15.7 mm |
|  | 0500-0468 | Sinker, Basket with Clip on Lid, 8 Mesh, 316 SS | 20.1 mm | 10.4 mm |
|  | 0500-0473 | Sinker, Helical with Epoxy Coating, 6 Turn Spiral, Magnetic | 21.3 mm | 9.4 mm |
|  | 0500-0474 | Sinker, 6 Turn Spiral, 316 SS | 21.3 mm | 9.4 mm |
|  | 0500-0475 | Sinker, XL, Teflon Coated | 29.2 mm | 11.8 mm |
|  | 0500-0476 | Sinker, "S" Shaped, 316 SS | 9.8 mm | 5.7 mm |
|  | 0500-0477 | Sinker, 5 Turn, 316 SS | 21.3 mm | 9.8 mm |
|  | 0500-0478 | Sinker, 4 Turn, SS, Inert Coated | 21.3 mm | 9.8 mm |
|  | 0500-0479 | Sinker, 5 Turn, SS, Inert Coated | 27.9 mm | 10.4 mm |
|  | 0500-0491 | Sinker, Swiss-Style, 316 SS | 15.5 mm | 5 mm |
|  | 0500-0492 | Sinker, Swiss-Style, 316 SS | 18 mm | 6 mm |
|  | 0500-0493 | Sinker, Swiss-Style, 316 SS | 19.3 mm | 7 mm |

CAPSULE SINKERS

| | PART NUMBER | DESCRIPTION | LENGTH | DIAMETER |
|---|-------------|---|---------|----------|
|  | 0500-0494 | Sinker, Swiss-Style, 316 SS | 23 mm | 8 mm |
|  | 0500-0495 | Sinker, Swiss-Style, 316 SS | 31 mm | 11 mm |
|  | 0500-0496 | Sinker, Japanese Pharmacopeia Basket, 316 SS | 23.7 mm | 11.9 mm |
|  | 5000-0011 | Sinker, 3 Prong Black Plastic | 25.4 mm | 5.1 mm |
|  | 5000-0012 | Sinker, 3 Prong Metallic Core, Inert Plastic Coated, Magnetic | 25.4 mm | 5.1 mm |

SINKER SIZE CHART

| PART NUMBER | SIZE 000 | SIZE 00 | SIZE 0 | SIZE AA | SIZE B | SUPRO A | SIZE 1 | SIZE 2 | SIZE 3 | SIZE 4 |
|-------------|----------|---------|--------|---------|--------|---------|--------|--------|--------|--------|
| 0500-0467 | | | | | | | | | | |
| 0500-0468 | | | | | | | | | | |
| 0500-0473 | | | | | | | | | | |
| 0500-0474 | | | | | | | | | | |
| 0500-0475 | | | | | | | | | | |
| 0500-0477 | | | | | | | | | | |
| 0500-0478 | | | | | | | | | | |
| 0500-0479 | | | | | | | | | | |

SHADED AREAS INDICATE THAT A CAPSULE FITS IN THE CORRESPONDING SINKER.



888.2.DISTEK
distekinc.com • info@distekinc.com
 121 North Center Drive • North Brunswick, NJ 08902