


Column Connection



Purchase the following connection kit when directly connecting the column to the injector of a nano-LC instrument.

MonoCap C18 HighResolution 2000 Connection Kit

Description	Contents	Cat.No.
MonoCap C18 HighResolution 2000 Connection Kit	1/16" PEEK Ferrule, SUS Nut, Sleeve, 2 pcs each. 1/32" PEEK Ferrule, SUS Nut, Sleeve, 2 pcs each.	5020-10017



Example of Column Connection




Column Connection

Purchase the following union when using a tubing that generates any indirect connection to the column.



Description	Model	Remarks	Cat.No.
Zero Dead Volume Union	U-435	For 1/16" SUS Tubing	6010-72352



Example of Column Connection using nanoViper



Example of Column Connection



Ordering Information

MonoCap™ C18 High Resolution 2000

I.D.(mm)	Length(mm)	Cat.No.
0.1	2000	5020-10015

- * A column stand is included.
- * End-fittings are not included.



MonoCap High Resolution 2000

Connection Kit for MonoCap™ C18 High Resolution



Connection Kit for MonoCap™ C18 High Resolution 2000

A dedicated connection kit for MonoCap C18 High Resolution 2000.

Use this connection kit when connecting the column directly to the system.

Description	Contents of Kit	Cat.No.
Connection Kit for MonoCap C18 High Resolution 2000	1/16" PEEK Ferrule, SUS Nut, Sleeve 2 pcs each 1/32" PEEK Ferrule, SUS Nut, Sleeve 2 pcs each	5020-10017

Zero Dead Volume Union

Connect the tubing from the system to this union and install the column to achieve zero dead volume.



Description	Model	Remarks	Cat.No.
Zero Dead Volume Union	U-435	For 1/16" SUS Tubing	6010-72352